

A20436

Leader in Biomolecular Solutions for Life Science



PAQR7 Rabbit pAb

Catalog No.: A20436

Basic Information

Observed MW

Refer to figures

Calculated MW

40kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to enable steroid binding activity and steroid hormone receptor activity. Predicted to be involved in response to steroid hormone. Predicted to be integral component of membrane. Predicted to be active in plasma membrane.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

164091

Swiss Prot

Q86WK9

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-60 of human PAQR7 (NP_848509.1).

Synonyms

mSR; MPRA; PGLP; PAQR7

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.