

A20437

Leader in Biomolecular Solutions for Life Science



# NREP Rabbit pAb

Catalog No.: A20437

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

8kDa

### Category

Polyclonal Antibody

### Applications

WB, ELISA

### Cross-Reactivity

Human

## Background

Predicted to be involved in axon regeneration; regulation of neuron differentiation; and regulation of transforming growth factor beta receptor signaling pathway. Predicted to be active in cytoplasm and nucleus.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

9315

### Swiss Prot

Q16612

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-68 of human NREP (NP\_004763.1).

### Synonyms

P311; PTZ17; SEZ17; D4S114; C5orf13; PRO1873; NREP

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.