

A20597

Leader in Biomolecular Solutions for Life Science



MRAP Rabbit pAb

Catalog No.: A20597

Basic Information

Observed MW

Refer to figures

Calculated MW

19kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene encodes a melanocortin receptor-interacting protein. The encoded protein regulates trafficking and function of the melanocortin 2 receptor in the adrenal gland. The encoded protein can also modulate signaling of other melanocortin receptors. Mutations in this gene have been associated with familial glucocorticoid deficiency type 2. Alternatively spliced transcript variants have been described.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

56246

Swiss Prot

Q8TCY5

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MRAP (NP_848932.1).

Synonyms

B27; FALP; FGD2; GCCD2; C21orf61; MRAP

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.