

A20600

Leader in Biomolecular Solutions for Life Science



## RFWD3 Rabbit pAb

Catalog No.: A20600

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

85kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human

### Background

Enables MDM2/MDM4 family protein binding activity; p53 binding activity; and ubiquitin protein ligase activity. Involved in several processes, including DNA metabolic process; regulation of cell cycle phase transition; and response to ionizing radiation. Located in nucleoplasm and site of DNA damage. Colocalizes with site of double-strand break. Implicated in Fanconi anemia.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

55159

#### Swiss Prot

Q6PCD5

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 150-450 of human RFWD3 (NP\_060594.3).

#### Synonyms

FANCW; RNF201; RFWD3

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.