

A20618

Leader in Biomolecular Solutions for Life Science



B3GAT3 Rabbit pAb

Catalog No.: A20618

1 Publications

Basic Information

Observed MW

Refer to figures

Calculated MW

37kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

The protein encoded by this gene is a member of the glucuronyltransferase gene family, enzymes that exhibit strict acceptor specificity, recognizing nonreducing terminal sugars and their anomeric linkages. This gene product catalyzes the formation of the glycosaminoglycan-protein linkage by way of a glucuronyl transfer reaction in the final step of the biosynthesis of the linkage region of proteoglycans. A pseudogene of this gene has been identified on chromosome 3.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

26229

Swiss Prot

O94766

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 29-190 of human B3GAT3 (NP_036332.2).

Synonyms

JDSCD; GLCATI; glcUAT-I; B3GAT3

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.