

A20645

Leader in Biomolecular Solutions for Life Science



## PDZK1IP1 Rabbit pAb

Catalog No.: A20645

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

12kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human

### Background

Located in extracellular exosome.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

10158

#### Swiss Prot

Q13113

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-114 of human PDZK1IP1 (NP\_005755.1).

#### Synonyms

DD96; SPAP; MAP17; PDZK1IP1

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.