

A20880

Leader in Biomolecular Solutions for Life Science



ZHX2 Rabbit mAb

Catalog No.: A20880

Recombinant

Basic Information

Observed MW

110kDa

Calculated MW

92kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

CloneNo number

ARC2826

Background

The members of the zinc fingers and homeoboxes gene family are nuclear homodimeric transcriptional repressors that interact with the A subunit of nuclear factor-Y (NF-YA) and contain two C2H2-type zinc fingers and five homeobox DNA-binding domains. This gene encodes member 2 of this gene family. In addition to forming homodimers, this protein heterodimerizes with member 1 of the zinc fingers and homeoboxes family.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

22882

Swiss Prot

Q9Y6X8

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 260-500 of human ZHX2 (Q9Y6X8).

Synonyms

RAF; AFR1; ZHX2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates, using ZHX2 antibody (A20880) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Enhanced Kit (RM00021).
Exposure time: 45s.