

A21085

Leader in Biomolecular Solutions for Life Science



Caspase-1 Rat mAb

Catalog No.: A21085

1 Publications

Basic Information

Observed MW

48kDa/33kDa

Calculated MW

46kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse, Rat

CloneNo number

ATC0001

Background

Enables cysteine-type endopeptidase activity. Involved in several processes, including positive regulation of I-kappaB kinase/NF-kappaB signaling; positive regulation of interleukin-1 beta production; and protein autoprocessing. Acts upstream of or within several processes, including membrane hyperpolarization; mitochondrial depolarization; and positive regulation of interleukin-1 alpha production. Located in extracellular region. Part of IPAF inflammasome complex. Is expressed in knee joint and metanephros. Orthologous to human CASP1 (caspase 1).

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

12362

Swiss Prot

P29452

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 206-220aa of mouse Caspase-1.

Synonyms

ICE; IL1bc

Contact



www.abclonal.com

Product Information

Source

Rat

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using Caspase-1 Rat mAb (A21085) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rat IgG (H+L) (AS028) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 180s.