

A21178

Leader in Biomolecular Solutions for Life Science



Siglec-8 Rabbit mAb

Catalog No.: A21178

Recombinant

Basic Information

Observed MW

Calculated MW

Category

SMab Recombinant Monoclonal Antibody

Applications

FC,ELISA

Cross-Reactivity

Human

CloneNo number

ARC52423

Background

Recommended Dilutions

FC 1:50 - 1:200

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 17-363 of human Siglec-8 (NP_055257.2).

Synonyms

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.