

A21224

Leader in Biomolecular Solutions for Life Science



CD79B Rabbit mAb

Catalog No.: A21224

Recombinant

Basic Information

Observed MW

30-50kDa

Calculated MW

26kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB, IHC-P, ELISA

Cross-Reactivity

Human, Mouse, Rat

CloneNo number

ARC52683

Background

The B lymphocyte antigen receptor is a multimeric complex that includes the antigen-specific component, surface immunoglobulin (Ig). Surface Ig non-covalently associates with two other proteins, Ig-alpha and Ig-beta, which are necessary for expression and function of the B-cell antigen receptor. This gene encodes the Ig-beta protein of the B-cell antigen component. Alternatively spliced transcript variants encoding different isoforms have been described.

Recommended Dilutions

WB	1:2000 - 1:20000
IHC-P	1:100 - 1:500

Immunogen Information

Gene ID

974

Swiss Prot

P40259

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 130-229 of human CD79B (NP_000617.1).

Synonyms

B29; IGB; AGM6; Igbeta; CD79B

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.


Validation Data




Western blot analysis of various lysates, using CD79B Rabbit mAb (A21224) at 1:20000 dilution.
 Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
 Lysates/proteins: 25µg per lane.
 Blocking buffer: 3% nonfat dry milk in TBST.
 Detection: ECL Basic Kit (RM00020).
 Exposure time: 180s.




Immunohistochemistry analysis of CD79B in paraffin-embedded Mouse spleen using CD79B Rabbit mAb (A21224) at dilution of 1:400(40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.




Immunohistochemistry analysis of CD79B in paraffin-embedded Human tonsil using CD79B Rabbit mAb (A21224) at dilution of 1:400(40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunohistochemistry analysis of CD79B in paraffin-embedded Rat spleen using CD79B Rabbit mAb (A21224) at dilution of 1:400(40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunohistochemistry analysis of CD79B in paraffin-embedded Human B-cell lymphoma using CD79B Rabbit mAb (A21224) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunohistochemistry analysis of CD79B in paraffin-embedded Human mantle cell lymphoma using CD79B Rabbit mAb (A21224) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.