

A21362

Leader in Biomolecular Solutions for Life Science



## CREB1 Rabbit pAb

Catalog No.: A21362

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

35kDa

#### Category

Polyclonal Antibody

#### Applications

IHC-P, IF/ICC, ELISA

#### Cross-Reactivity

Human, Mouse, Rat

### Background

This gene encodes a transcription factor that is a member of the leucine zipper family of DNA binding proteins. This protein binds as a homodimer to the cAMP-responsive element, an octameric palindrome. The protein is phosphorylated by several protein kinases, and induces transcription of genes in response to hormonal stimulation of the cAMP pathway. Alternate splicing of this gene results in several transcript variants encoding different isoforms.

### Recommended Dilutions

IHC-P	1:50 - 1:200
IF/ICC	1:50 - 1:200

### Immunogen Information

Gene ID	Swiss Prot
1385	P16220

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-341 of human CREB1 (NP\_604391.1).

#### Synonyms

CREB; CREB-1; CREB1

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information


Source	Isotype	Purification
Rabbit	IgG	Affinity purification

#### Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.


## Validation Data




Immunohistochemistry analysis of CREB1 in paraffin-embedded human colon carcinoma using CREB1 Rabbit pAb (A21362) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.




Immunohistochemistry analysis of CREB1 in paraffin-embedded human esophageal cancer using CREB1 Rabbit pAb (A21362) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.




Immunohistochemistry analysis of CREB1 in paraffin-embedded mouse brain using CREB1 Rabbit pAb (A21362) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunohistochemistry analysis of CREB1 in paraffin-embedded rat brain using CREB1 Rabbit pAb (A21362) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunohistochemistry analysis of CREB1 in paraffin-embedded rat ovary using CREB1 Rabbit pAb (A21362) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunofluorescence analysis of PC-12 cells using CREB1 Rabbit pAb (A21362) at dilution of 1:50 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining.