

A21431

Leader in Biomolecular Solutions for Life Science



## DNAJC2 Rabbit pAb

Catalog No.: A21431

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

72kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human, Mouse

### Background

This gene is a member of the M-phase phosphoprotein (MPP) family. The gene encodes a phosphoprotein with a J domain and a Myb DNA-binding domain which localizes to both the nucleus and the cytosol. The protein is capable of forming a heterodimeric complex that associates with ribosomes, acting as a molecular chaperone for nascent polypeptide chains as they exit the ribosome. This protein was identified as a leukemia-associated antigen and expression of the gene is upregulated in leukemic blasts. Also, chromosomal aberrations involving this gene are associated with primary head and neck squamous cell tumors. This gene has a pseudogene on chromosome 6. Alternatively spliced variants which encode different protein isoforms have been described.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

27000

#### Swiss Prot

Q99543

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human DNAJC2 (NP\_055192.1).

#### Synonyms

ZRF1; ZUO1; MPP11; MPHOSPH11; DNAJC2

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification


Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using DNAJC2 Rabbit pAb (A21431) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 30s.