

A21855

Leader in Biomolecular Solutions for Life Science



[KO Validated] PDCD4 Rabbit mAb

Catalog No.: A21855 **KO** **Validated** **Recombinant**

Basic Information

Observed MW

Refer to figures

Calculated MW

52kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene is a tumor suppressor and encodes a protein that binds to the eukaryotic translation initiation factor 4A1 and inhibits its function by preventing RNA binding. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

27250

Swiss Prot

Q53EL6

Immunogen

Recombinant protein of human PDCD4.

Synonyms

H731; D4

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.