

A21857

Leader in Biomolecular Solutions for Life Science



## [KO Validated] CD168/RHAMM Rabbit mAb

Catalog No.: A21857

**KO** Validated

Recombinant

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

84kDa

#### Category

SMab Recombinant Monoclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human

### Background

The protein encoded by this gene is involved in cell motility. It is expressed in breast tissue and together with other proteins, it forms a complex with BRCA1 and BRCA2, thus is potentially associated with higher risk of breast cancer. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

3161

#### Swiss Prot

O75330

#### Immunogen

Recombinant protein of human CD168/RHAMM.

#### Synonyms

CD168; IHABP; RHAMM; MM

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.