

A21859

Leader in Biomolecular Solutions for Life Science



[KO Validated] LAMP2 Rabbit mAb

Catalog No.: A21859

KO Validated

Recombinant

Basic Information

Observed MW

Refer to figures

Calculated MW

45kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may play a role in tumor cell metastasis. It may also function in the protection, maintenance, and adhesion of the lysosome. Alternative splicing of this gene results in multiple transcript variants encoding distinct proteins.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

3920

Swiss Prot

P13473

Immunogen

Recombinant protein of human LAMP2.

Synonyms

DND; LAMPB; CD107b; LAMP-2; LGP-96; LGP110; P2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.