

A21863

Leader in Biomolecular Solutions for Life Science



[KO Validated] NBS1/NBN Rabbit mAb

Catalog No.: A21863

KO Validated

Recombinant

Basic Information

Observed MW

Refer to figures

Calculated MW

85kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Mutations in this gene are associated with Nijmegen breakage syndrome, an autosomal recessive chromosomal instability syndrome characterized by microcephaly, growth retardation, immunodeficiency, and cancer predisposition. The encoded protein is a member of the MRE11/RAD50 double-strand break repair complex which consists of 5 proteins. This gene product is thought to be involved in DNA double-strand break repair and DNA damage-induced checkpoint activation.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

4683

Swiss Prot

O60934

Immunogen

Recombinant protein of human NBS1/NBN.

Synonyms

ATV; NBS; P95; NBS1; AT-V1; AT-V2; BN

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.