

A21867

Leader in Biomolecular Solutions for Life Science



## [KO Validated] NSD3/WHSC1L1 Rabbit mAb

Catalog No.: A21867

**KO** Validated

Recombinant

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

162kDa

#### Category

SMab Recombinant Monoclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human

### Background

This gene is related to the Wolf-Hirschhorn syndrome candidate-1 gene and encodes a protein with PWWP (proline-tryptophan-tryptophan-proline) domains. This protein methylates histone H3 at lysine residues 4 and 27, which represses gene transcription. Two alternatively spliced variants have been described.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

54904

#### Swiss Prot

Q9BZ95

#### Immunogen

Recombinant protein of human NSD3/WHSC1L1.

#### Synonyms

KMT3F; KMT3G; WHISTLE; WHSC1L1; pp14328; L1

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.