

A21883

Leader in Biomolecular Solutions for Life Science



## [KO Validated] GNAI2 Rabbit mAb

Catalog No.: A21883

**KO** Validated

Recombinant

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

40kDa

#### Category

SMab Recombinant Monoclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human

### Background

The protein encoded by this gene is an alpha subunit of guanine nucleotide binding proteins (G proteins). The encoded protein contains the guanine nucleotide binding site and is involved in the hormonal regulation of adenylate cyclase. Several transcript variants encoding different isoforms have been found for this gene.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

2771

#### Swiss Prot

P04899

#### Immunogen

Recombinant protein of human GNAI2.

#### Synonyms

GIP; HG1C; GNAI2B; H\_LUCA15.1; H\_LUCA16.1; I2

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.