

A22598

Leader in Biomolecular Solutions for Life Science



α -Synuclein Rabbit PolymAb®

Catalog No.: A22598

Basic Information

Observed MW

18kDa

Calculated MW

14kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB, IHC-P, IF/ICC, ELISA

Cross-Reactivity

Human, Mouse, Rat

CloneNo number

ARC50333_ARC50335

Background

Alpha-synuclein is a member of the synuclein family, which also includes beta- and gamma-synuclein. Synucleins are abundantly expressed in the brain and alpha- and beta-synuclein inhibit phospholipase D2 selectively. SNCA may serve to integrate presynaptic signaling and membrane trafficking. Defects in SNCA have been implicated in the pathogenesis of Parkinson disease. SNCA peptides are a major component of amyloid plaques in the brains of patients with Alzheimer's disease. Alternatively spliced transcripts encoding different isoforms have been identified for this gene.

Recommended Dilutions

WB	1:500 - 1:1000
IHC-P	1:1000-1:5000
IF/ICC	1:50 - 1:200

Immunogen Information

Gene ID

6622

Swiss Prot

P37840

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 59-140 of human α -Synuclein Rabbit PolymAb® (NP_000336.1).

Synonyms

PD1; NACP; PARK1; PARK4

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.



Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.

Validation Data




Western blot analysis of various lysates, using α -Synuclein Rabbit PolymAb® (A22598) at 1:1000 dilution.
 Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
 Lysates/proteins: 25 μ g per lane.
 Blocking buffer: 3% nonfat dry milk in TBST.
 Detection: ECL Basic Kit (RM00020).
 Exposure time: 10s.


Western blot analysis of various lysates, using α -Synuclein Rabbit PolymAb® (A22598) at 1:1000 dilution.
 Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
 Lysates/proteins: 25 μ g per lane.
 Blocking buffer: 3% nonfat dry milk in TBST.
 Detection: ECL Basic Kit (RM00020).
 Exposure time: 90s.



Immunohistochemistry analysis of α -Synuclein in paraffin-embedded mouse brain tissue using α -Synuclein Rabbit PolymAb® (A22598) at a dilution of 1:4000 (40x lens). High pressure antigen retrieval was performed with 0.01 M citrate buffer (pH 6.0) prior to IHC staining.




Immunohistochemistry analysis of α -Synuclein in paraffin-embedded human brain tissue using α -Synuclein Rabbit PolymAb® (A22598) at a dilution of 1:4000 (40x lens). High pressure antigen retrieval was performed with 0.01 M citrate buffer (pH 6.0) prior to IHC staining.



Immunohistochemistry analysis of α -Synuclein in paraffin-embedded rat brain tissue using α -Synuclein Rabbit PolymAb® (A22598) at a dilution of 1:4000 (40x lens). High pressure antigen retrieval was performed with 0.01 M citrate buffer (pH 6.0) prior to IHC staining.

Validation Data



Immunofluorescence analysis of Neuro-2a cells using α -Synuclein Rabbit PolymAb[®] (A22598) at dilution of 1:50 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining.