

A23075

Leader in Biomolecular Solutions for Life Science



# EpCAM Rabbit mAb

Catalog No.: A23075

Recombinant

## Basic Information

### Observed MW

40kDa

### Calculated MW

35kDa

### Category

SMab Recombinant Monoclonal Antibody

### Applications

WB, IHC-P, ELISA

### Cross-Reactivity

Human

### CloneNo number

ARC57675

## Background

This gene encodes a carcinoma-associated antigen and is a member of a family that includes at least two type I membrane proteins. This antigen is expressed on most normal epithelial cells and gastrointestinal carcinomas and functions as a homotypic calcium-independent cell adhesion molecule. The antigen is being used as a target for immunotherapy treatment of human carcinomas. Mutations in this gene result in congenital tufting enteropathy.

## Recommended Dilutions

WB 1:10000 - 1:40000

IHC-P 1:50 - 1:200

## Immunogen Information

### Gene ID

4072

### Swiss Prot

P16422

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 24-265 of human EpCAM (NP\_002345.2).

### Synonyms

ESA; KSA; M4S1; MK-1; DIAR5; EGP-2; EGP40; KS1/4; MIC18; TROP1; BerEp4; EGP314; HNPCC8; LYNCH8; MOC-31; Ber-Ep4; TACSTD1; EpCAM

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification


Affinity purification

### Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.


## Validation Data



Western blot analysis of various lysates, using EpCAM Rabbit mAb (A23075) at 1:40000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 180s.



Immunohistochemistry analysis of EpCAM in paraffin-embedded human colon carcinoma using EpCAM Rabbit mAb (A23075) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunohistochemistry analysis of EpCAM in paraffin-embedded human kidney using EpCAM Rabbit mAb (A23075) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.