

A23592

Leader in Biomolecular Solutions for Life Science



PE Rabbit anti-Dog CD19 mAb

Catalog No.: A23592

Basic Information

Observed MW

Calculated MW

Category

SMab Recombinant Monoclonal Antibody

Applications

FC

Cross-Reactivity

Dog

CloneNo number

ARC58529-PE

Conjugate

PE. Ex:565nm. Em:574nm.

Background

Recommended Dilutions

FC 5 µl per 10⁶ cells in
 100 µl volume

Immunogen Information

Gene ID

607898

Swiss Prot

Immunogen

Recombinant protein of dog CD19.

Synonyms

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification



Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.03% proclin300,0.2% BSA,pH7.3.

Validation Data



Flow cytometry: 1×10^6 293F cells(negative control,left) and 293F(Transfection,right) cells were surface-stained with PE Rabbit anti-Dog CD19 mAb(A23592,5 μ l/Test,orange line) or PE Rabbit IgG isotype control (5 μ l/Test,blue line).Non-fluorescently stained cells was used as blank control (red line).

Flow cytometry: 1×10^6 293F(Transfection) cells were surface-stained with PE Rabbit IgG isotype control (5 μ l/Test,left) or PE Rabbit anti-Dog CD19 mAb(A23592,5 μ l/Test,right).