

A23819

Leader in Biomolecular Solutions for Life Science



## ABflo® 647 Rabbit anti-Pig CD3G mAb

Catalog No.: A23819

### Basic Information

**Observed MW**

**Calculated MW**

**Category**

SMab Recombinant Monoclonal Antibody

**Applications**

FC

**Cross-Reactivity**

Pig

**CloneNo number**

ARC60516-ABflo647

**Conjugate**

ABflo® 647. Ex:648nm. Em:664nm.

### Background

### Recommended Dilutions

**FC**                      5 µl per 10<sup>6</sup> cells in  
                                 100 µl volume

### Immunogen Information

**Gene ID**

494013

**Swiss Prot**

**Immunogen**

Recombinant protein of pig CD3G.

**Synonyms**

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**


Affinity purification

**Storage**


Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.03% proclin300,0.2% BSA,pH7.3.

## Validation Data



Flow cytometry:  $1 \times 10^6$  293F cells (negative control, left) and 293F (Transfection, right) cells were surface-stained with ABflo® 647 Rabbit anti-Pig CD3G mAb (A23819, 5  $\mu$ l/Test, orange line) or ABflo® 647 Rabbit IgG isotype control (A22070, 5  $\mu$ l/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).



Flow cytometry:  $1 \times 10^6$  293F Pig CD3G cells were surface-stained with ABflo® 647 Rabbit IgG isotype control (A22070, 5  $\mu$ l/Test, left) or ABflo® 647 Rabbit anti-Pig CD3G mAb (A23819, 5  $\mu$ l/Test, right).