

Pan Acetyl-Lysine pAb

Catalog No.	A2391	Category	Acetylated Antibodies
Applications	WB, IP	Observed MW	Refer to Figures
Cross-reactivity	ALL	Calculated MW	

Immunogen Information

Immunogen Acetylated BSA

Gene ID

Swiss prot

Synonyms

Product information

Source Rabbit

Isotype IgG

Purification method Affinity purification


Storage Store at -20°C. Avoid freeze / thaw cycles.
 Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 - 1:1000

IP 1:20 - 1:50



Western blot analysis of extracts of HeLa cells, using Pan Acetyl-Lysine antibody (A2391).

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.
 Blocking buffer: 3% nonfat dry milk in TBST.