

A24261

Leader in Biomolecular Solutions for Life Science



PP2A-B55 α /PR55 α /PPP2R2A Rabbit mAb

Catalog No.: A24261

Recombinant

Basic Information

Observed MW

52kDa

Calculated MW

52kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

CloneNo number

ARC60262

Background

The product of this gene belongs to the phosphatase 2 regulatory subunit B family. Protein phosphatase 2 is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. The B regulatory subunit might modulate substrate selectivity and catalytic activity. This gene encodes an alpha isoform of the regulatory subunit B55 subfamily. Alternatively spliced transcript variants have been described.

Recommended Dilutions

WB 1:1000 - 1:5000

Immunogen Information

Gene ID

5520

Swiss Prot

P63151

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 94-194 of human PP2A-B55 α /PR55 α /PPP2R2A (NP_002708.1).

Synonyms

B55A; PR52A; PR55A; B55ALPHA; PR55alpha; PP2A-B55 α /PR55 α /PPP2R2A

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of lysates from HeLa cells, using PP2A-B55α/PR55α/PPP2R2A Rabbit mAb (A24261) at 1:2000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25ug per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 30s.