

A24481

Leader in Biomolecular Solutions for Life Science



CK14 Rabbit pAb

Catalog No.: A24481

Basic Information

Observed MW

Refer to figures

Calculated MW

51kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene encodes a member of the keratin family, the most diverse group of intermediate filaments. This gene product, a type I keratin, is usually found as a heterotetramer with two keratin 5 molecules, a type II keratin. Together they form the cytoskeleton of epithelial cells. Mutations in the genes for these keratins are associated with epidermolysis bullosa simplex. At least one pseudogene has been identified at 17p12-p11.

Recommended Dilutions

WB 1:1000 - 1:5000

Immunogen Information

Gene ID

3861

Swiss Prot

P02533

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 372-472 of human CK14 (NP_000517.3).

Synonyms

KRT14; CK14; EBS3; EBS4; K14; NFJ; keratin 14

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates, using CK14 Rabbit pAb (A24481) at 1:3000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25ug per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Negative control (NC): SH-SY5Y
Exposure time: 5s.