

A2926

Leader in Biomolecular Solutions for Life Science



# FOXD3 Rabbit pAb

Catalog No.: A2926

## Basic Information

### Observed MW

48kDa

### Calculated MW

48kDa

### Category

Polyclonal Antibody

### Applications

WB

### Cross-Reactivity

Human, Mouse

## Background

This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. Mutations in this gene cause autoimmune susceptibility 1.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

27022

### Swiss Prot

Q9UJU5

### Immunogen

A synthetic peptide corresponding to a sequence within amino acids 322-335AA of human FOXD3.

### Synonyms

AIS1; HFH2; VAMAS2; Genesis; FOXD3

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Validation Data

---



Western blot analysis of lysates from HeLa cells, using FOXD3 Rabbit pAb (A2926).  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.