

A3052

Leader in Biomolecular Solutions for Life Science



NCAN Rabbit pAb

Catalog No.: A3052

Basic Information

Observed MW

Refer to figures

Calculated MW

143kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human, Mouse, Rat

Background

Neurocan is a chondroitin sulfate proteoglycan thought to be involved in the modulation of cell adhesion and migration.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

1463

Swiss Prot

O14594

Immunogen

A synthetic peptide of human NCAN

Synonyms

CSPG3; NCAN

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.