

A3705

Leader in Biomolecular Solutions for Life Science



KMT1B / SUV39H2 Rabbit mAb

Catalog No.: A3705

Recombinant

Basic Information

Observed MW

45kDa

Calculated MW

47kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

CloneNo number

ARC0829

Background

Enables S-adenosyl-L-methionine binding activity; histone methyltransferase activity (H3-K9 specific); and zinc ion binding activity. Involved in chromatin assembly or disassembly and chromatin remodeling. Acts upstream of or within cellular response to hypoxia and negative regulation of transcription by RNA polymerase II. Located in chromatin.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

79723

Swiss Prot

Q9H5I1

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 311-410 of human KMT1B / SUV39H2 (Q9H5I1).

Synonyms

KMT1B; KMT1B / SUV39H2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using KMT1B / SUV39H2 antibody (A3705) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 30s.