

A4370

Leader in Biomolecular Solutions for Life Science



EDIL3 Rabbit pAb

Catalog No.: A4370

Basic Information

Observed MW

54kDa

Calculated MW

54kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

The protein encoded by this gene is an integrin ligand. It plays an important role in mediating angiogenesis and may be important in vessel wall remodeling and development. It also influences endothelial cell behavior.

Recommended Dilutions

WB 1:1000 - 1:2000

Immunogen Information

Gene ID

10085

Swiss Prot

O43854

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 60-300 of human EDIL3 (NP_005702.3).

Synonyms

DEL1; EDIL3

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using EDIL3 antibody (A4370) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 1s.