

A4442

Leader in Biomolecular Solutions for Life Science



IVNS1ABP Rabbit pAb

Catalog No.: A4442

Basic Information

Observed MW

72kDa

Calculated MW

72kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Involved in RNA splicing; negative regulation of protein ubiquitination; and response to virus. Located in cytosol. Implicated in immunodeficiency 70.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

10625

Swiss Prot

Q9Y6Y0

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human IVNS1ABP (NP_006460.2).

Synonyms

ND1; ARA3; NS-1; IMD70; NS1BP; FLARA3; KLHL39; NS1-BP; HSPC068; IVNS1ABP

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using IVNS1ABP antibody (A4442) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 30s.