

A4555

Leader in Biomolecular Solutions for Life Science



Olig2 Rabbit mAb

Catalog No.: A4555

Recombinant

2 Publications

Basic Information

Observed MW

40kDa

Calculated MW

32kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse, Rat

CloneNo number

ARC1032

Background

This gene encodes a basic helix-loop-helix transcription factor which is expressed in oligodendroglial tumors of the brain. The protein is an essential regulator of ventral neuroectodermal progenitor cell fate. The gene is involved in a chromosomal translocation t(14;21)(q11.2;q22) associated with T-cell acute lymphoblastic leukemia. Its chromosomal location is within a region of chromosome 21 which has been suggested to play a role in learning deficits associated with Down syndrome.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

10215

Swiss Prot

Q13516

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Olig2 (Q13516).

Synonyms

BHLHB1; OLIGO2; RACK17; PRKCBP2; bHLHe19; Olig2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using Olig2 Rabbit mAb (A4555) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 3min.