

A5508

Leader in Biomolecular Solutions for Life Science



## TRIM22 Rabbit pAb

Catalog No.: A5508

### Basic Information

**Observed MW**

57kDa

**Calculated MW**

57kDa

**Category**

Polyclonal Antibody

**Applications**

ELISA

**Cross-Reactivity**

Human

### Background

The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. This protein localizes to the cytoplasm and its expression is induced by interferon. The protein is involved in innate immunity against different DNA and RNA viruses.

### Recommended Dilutions

### Immunogen Information

**Gene ID**

10346

**Swiss Prot**

Q8IYM9

**Immunogen**

Recombinant fusion protein containing a sequence corresponding to amino acids 215-411 of human TRIM22 (NP\_006065.2).

**Synonyms**

RNF94; STAF50; GPSTAF50; TRIM22

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**

Affinity purification

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data

