

A6325

Leader in Biomolecular Solutions for Life Science



CFDP1 Rabbit pAb

Catalog No.: A6325

Basic Information

Observed MW

33kDa

Calculated MW

34kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Predicted to act upstream of or within several processes, including cell adhesion; negative regulation of fibroblast apoptotic process; and regulation of cell shape. Predicted to be located in kinetochore.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

10428

Swiss Prot

Q9UEE9

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-299 of human CFDP1 (NP_006315.1).

Synonyms

p97; BCNT; CP27; SWC5; Yeti; CENP-29; BUCENTAUR; CFDP1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using CFDP1 antibody (A6325) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25 μ g per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Enhanced Kit (RM00021).
Exposure time: 90s.