

A6442

Leader in Biomolecular Solutions for Life Science



## UBE3C Rabbit pAb

Catalog No.: A6442

2 Publications

### Basic Information

#### Observed MW

124kDa

#### Calculated MW

124kDa

#### Category

Polyclonal Antibody

#### Applications

WB,IF/ICC,ELISA

#### Cross-Reactivity

Human

### Background

Enables ubiquitin protein ligase activity. Involved in protein polyubiquitination. Predicted to be located in nucleus. Predicted to be part of proteasome complex.

### Recommended Dilutions

WB	1:500 - 1:2000
IF/ICC	1:10 - 1:100

### Immunogen Information

#### Gene ID

9690

#### Swiss Prot

Q15386

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 970-1083 of human UBE3C (NP\_055486.2).

#### Synonyms

RAUL; HECTH2; NDSMBA; NEDSMBA; UBE3C

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.