

A6513

Leader in Biomolecular Solutions for Life Science



FOXF1 Rabbit mAb

Catalog No.: A6513

Recombinant

Basic Information

Observed MW

40kDa

Calculated MW

40kDa

Category

SMab Recombinant Monoclonal
Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse

CloneNo number

ARC1863

Background

This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in the regulation of pulmonary genes as well as embryonic development.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

2294

Swiss Prot

Q12946

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOXF1 (Q12946).

Synonyms

FKHL5; ACDMPV; FREAC1; FOXF1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using FOXF1 antibody (A6513) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 10s.