

A6590

Leader in Biomolecular Solutions for Life Science



FLOT2 Rabbit pAb

Catalog No.: A6590

1 Publications

Basic Information

Observed MW

47kDa

Calculated MW

47kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Caveolae are small domains on the inner cell membrane involved in vesicular trafficking and signal transduction. This gene encodes a caveolae-associated, integral membrane protein, which is thought to function in neuronal signaling.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

2319

Swiss Prot

Q14254

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 239-428 of human FLOT2 (NP_004466.2).

Synonyms

ESA; ECS1; ESA1; ECS-1; M17S1; FLOT2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using FLOT2 Rabbit pAb (A6590) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 90s.