

A6628

Leader in Biomolecular Solutions for Life Science



IL13RA1 Rabbit pAb

Catalog No.: A6628

1 Publications

Basic Information

Observed MW

53kDa

Calculated MW

49kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

The protein encoded by this gene is a subunit of the interleukin 13 receptor. This subunit forms a receptor complex with IL4 receptor alpha, a subunit shared by IL13 and IL4 receptors. This subunit serves as a primary IL13-binding subunit of the IL13 receptor, and may also be a component of IL4 receptors. This protein has been shown to bind tyrosine kinase TYK2, and thus may mediate the signaling processes that lead to the activation of JAK1, STAT3 and STAT6 induced by IL13 and IL4.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

3597

Swiss Prot

P78552

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 29-343 of human IL13RA1 (NP_001551.1).

Synonyms

NR4; CT19; CD213A1; IL-13Ra; IL13RA1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using IL13RA1 antibody (A6628) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 90s.