

A6885

Leader in Biomolecular Solutions for Life Science



## CKS1B Rabbit pAb

Catalog No.: A6885

### Basic Information

**Observed MW**

55kDa

**Calculated MW**

10kDa

**Category**

Polyclonal Antibody

**Applications**

WB, ELISA

**Cross-Reactivity**

Human

### Background

CKS1B protein binds to the catalytic subunit of the cyclin dependent kinases and is essential for their biological function. The CKS1B mRNA is found to be expressed in different patterns through the cell cycle in HeLa cells, which reflects a specialized role for the encoded protein. At least two transcript variants have been identified for this gene, and it appears that only one of them encodes a protein.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

**Gene ID**

1163

**Swiss Prot**

P61024

**Immunogen**

Recombinant fusion protein containing a sequence corresponding to amino acids 1-79 of human CKS1B (NP\_001817.1).

**Synonyms**

CKS1; ckshs1; PNAS-16; PNAS-18; CKS1B

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**

Affinity purification

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.