

A7053

Leader in Biomolecular Solutions for Life Science



GJC1 Rabbit pAb

Catalog No.: A7053

1 Publications

Basic Information

Observed MW

45kDa

Calculated MW

45kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

This gene is a member of the connexin gene family. The encoded protein is a component of gap junctions, which are composed of arrays of intercellular channels that provide a route for the diffusion of low molecular weight materials from cell to cell. Alternatively spliced transcript variants encoding the same isoform have been described.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

10052

Swiss Prot

P36383

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 247-396 of human GJC1 (NP_005488.2).

Synonyms

CX45; GJA7; GJC1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using GJC1 antibody (A7053) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 30s.