

A7151

Leader in Biomolecular Solutions for Life Science



## U2AF1L4 Rabbit pAb

Catalog No.: A7151

### Basic Information

**Observed MW**

Refer to figures

**Calculated MW**

26kDa

**Category**

Polyclonal Antibody

**Applications**

WB, ELISA

**Cross-Reactivity**

Human

### Background

Predicted to enable pre-mRNA 3'-splice site binding activity. Predicted to be involved in mRNA splicing, via spliceosome. Predicted to be located in nucleoplasm. Predicted to be part of U2AF complex and spliceosomal complex.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

**Gene ID**

199746

**Swiss Prot**

Q8WU68

**Immunogen**

Recombinant fusion protein containing a sequence corresponding to amino acids 1-202 of human U2AF1L4 (NP\_659424.2).

**Synonyms**

U2af26; U2AF1L3; U2AF1RS3; U2AF1-RS3; U2AF1L3V1; U2AF1L4

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**

Affinity purification

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.