

A7268

Leader in Biomolecular Solutions for Life Science



LMO4 Rabbit pAb

Catalog No.: A7268

Basic Information

Observed MW

Calculated MW

18kDa

Category

Polyclonal Antibody

Applications

IF/ICC,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

This gene encodes a cysteine-rich protein that contains two LIM domains but lacks a DNA-binding homeodomain. The encoded protein may play a role as a transcriptional regulator or as an oncogene.

Recommended Dilutions

IF/ICC 1:50 - 1:100

Immunogen Information

Gene ID

8543

Swiss Prot

P61968

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human LMO4 (NP_006760.1).

Synonyms

LMO4

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Immunofluorescence analysis of MCF-7 cells using LMO4 antibody (A7268). Blue: DAPI for nuclear staining.