

A7478

Leader in Biomolecular Solutions for Life Science



## ZFAND3 Rabbit pAb

Catalog No.: A7478

2 Publications

### Basic Information

**Observed MW**

27kDa

**Calculated MW**

25kDa

**Category**

Polyclonal Antibody

**Applications**

IHC-P,ELISA

**Cross-Reactivity**

Human,Mouse,Rat

### Background

Predicted to enable DNA binding activity and zinc ion binding activity.

### Recommended Dilutions

IHC-P 1:50 - 1:200

### Immunogen Information

**Gene ID**

60685

**Swiss Prot**

Q9H8U3

**Immunogen**

Recombinant fusion protein containing a sequence corresponding to amino acids 1-227 of human ZFAND3 (NP\_068762.1).

**Synonyms**

TEX27; ZFAND3

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**

Affinity purification


**Storage**

Store at -20°C. Avoid freeze / thaw cycles.


Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Validation Data

---



Immunohistochemistry analysis of paraffin-embedded rat kidney using ZFAND3 antibody (A7478) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.



Immunohistochemistry analysis of paraffin-embedded mouse kidney using ZFAND3 antibody (A7478) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.