

A7506

Leader in Biomolecular Solutions for Life Science



RAB31 Rabbit pAb

Catalog No.: A7506

1 Publications

Basic Information

Observed MW

22kDa

Calculated MW

22kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Enables GDP binding activity and GTP binding activity. Involved in several processes, including Golgi to plasma membrane protein transport; cellular response to insulin stimulus; and receptor internalization. Located in early endosome; phagocytic vesicle; and trans-Golgi network membrane. Biomarker of severe acute respiratory syndrome.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

11031

Swiss Prot

Q13636

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-195 of human RAB31 (NP_006859.2).

Synonyms

Rab22B; RAB31

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using RAB31 antibody (A7506) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 30s.