

A7583

Leader in Biomolecular Solutions for Life Science



WDR47 Rabbit pAb

Catalog No.: A7583

Basic Information

Observed MW

Refer to figures

Calculated MW

102kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Predicted to be located in cytoplasm and microtubule.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

22911

Swiss Prot

O94967

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 350-600 of human WDR47 (NP_001136023.1).

Synonyms

WDR47

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.