

A7903

Leader in Biomolecular Solutions for Life Science



IFNA6 Rabbit pAb

Catalog No.: A7903

Basic Information

Observed MW

Refer to figures

Calculated MW

22kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to enable cytokine activity and type I interferon receptor binding activity. Predicted to be involved in several processes, including B cell activation; lymphocyte activation involved in immune response; and positive regulation of peptidyl-serine phosphorylation of STAT protein. Predicted to be located in extracellular region. Predicted to be active in extracellular space. Biomarker of anogenital venereal wart.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

3443

Swiss Prot

P05013

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 21-189 of human IFNA6 (NP_066282.1).

Synonyms

IFN-alphaK; IFNA6

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.