

A8124

Leader in Biomolecular Solutions for Life Science



IFNA4 Rabbit pAb

Catalog No.: A8124

Basic Information

Observed MW

Refer to figures

Calculated MW

22kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to enable cytokine activity and type I interferon receptor binding activity. Predicted to be involved in several processes, including B cell activation; lymphocyte activation involved in immune response; and positive regulation of peptidyl-serine phosphorylation of STAT protein. Predicted to be located in extracellular region. Predicted to be active in extracellular space.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

3441

Swiss Prot

P05014

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 24-189 of human IFNA4 (NP_066546.1).

Synonyms

INFA4; IFN-alpha4a; IFNA4

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.