

A8336

Leader in Biomolecular Solutions for Life Science



DROSHA Rabbit pAb

Catalog No.: A8336

Basic Information

Observed MW

159kDa

Calculated MW

159kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

This gene encodes a ribonuclease (RNase) III double-stranded RNA-specific ribonuclease and subunit of the microprocessor protein complex, which catalyzes the initial processing step of microRNA (miRNA) synthesis. The encoded protein cleaves the stem loop structure from the primary microRNA (pri-miRNA) in the nucleus, yielding the precursor miRNA (pre-miRNA), which is then exported to the cytoplasm for further processing. In a human cell line lacking a functional copy of this gene, canonical miRNA synthesis is reduced. Somatic mutations in this gene have been observed in human patients with kidney cancer.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

29102

Swiss Prot

Q9NRR4

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DROSHA (NP_037367.3).

Synonyms

RN3; ETOH12; RNASEN; RANSE3L; RNASE3L; HSA242976; DROSHA

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of lysates from Mouse heart, using DROSHA Rabbit pAb (A8336) at 1:500 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 180s.