

A8954

Leader in Biomolecular Solutions for Life Science



## C1QTNF4 Rabbit pAb

Catalog No.: A8954

### Basic Information

**Observed MW**

38kDa

**Calculated MW**

35kDa

**Category**

Polyclonal Antibody

**Applications**

WB, ELISA

**Cross-Reactivity**

Human

### Background

Predicted to enable cytokine activity. Involved in positive regulation of cytokine production and positive regulation of signal transduction. Located in extracellular space.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

**Gene ID**

114900

**Swiss Prot**

Q9BXJ3

**Immunogen**

Recombinant fusion protein containing a sequence corresponding to amino acids 250-329 of human C1QTNF4 (NP\_114115.2).

**Synonyms**

CTRP4; ZACRP4; C1QTNF4

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**


Affinity purification

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using C1QTNF4 Rabbit pAb (A8954) at 1:3000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 5s.