

A9165

Leader in Biomolecular Solutions for Life Science



Fyn Rabbit mAb

Catalog No.: A9165

Recombinant

1 Publications

Basic Information

Observed MW

59kDa

Calculated MW

61kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

CloneNo number

ARC1459

Background

This gene is a member of the protein-tyrosine kinase oncogene family. It encodes a membrane-associated tyrosine kinase that has been implicated in the control of cell growth. The protein associates with the p85 subunit of phosphatidylinositol 3-kinase and interacts with the fyn-binding protein. Alternatively spliced transcript variants encoding distinct isoforms exist.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

2534

Swiss Prot

P06241

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Fyn (P06241).

Synonyms

SLK; SYN; p59-FYN; Fyn

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using Fyn Rabbit mAb (A9165) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 1min.