

A9338

Leader in Biomolecular Solutions for Life Science



FUT7 Rabbit pAb

Catalog No.: A9338

Basic Information

Observed MW

50kDa

Calculated MW

39kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

The protein encoded by this gene is a Golgi stack membrane protein that is involved in the creation of sialyl-Lewis X antigens. The encoded protein can direct the synthesis of the E-selectin-binding sialyl-Lewis X moiety.

Recommended Dilutions

WB 1:200 - 1:2000

Immunogen Information

Gene ID

2529

Swiss Prot

Q11130

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 100-200 of human FUT7 (NP_004470.1).

Synonyms

FucT-VII; FUT7

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using FUT7 Antibody (A9338) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 15s.